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Fast and Precise Isotope Analysis of Liquids on a DELTA V Isotope Ratio MS with High Temperature Conversion Elemental Analyzer

Oliver Kracht, Andreas Hilkert, Thermo Fisher Scientific, Bremen, Germany

Key Words

- TC/EA
- Isotope Ratio MS
- EA-IRMS
- Water Analysis

Introduction

^{18}O and ^2H isotope ratio analysis of aqueous solutions such as water, urine, blood plasma, wine, can be performed by a variety of different techniques: equilibration, chromium reduction and carbon reduction (also referred to as high temperature conversion or pyrolysis). The high temperature conversion (TC) method allows a direct isotope ratio analysis of both oxygen and hydrogen in continuous flow mode. It offers the analysis of two isotopes in five minutes from only one injection of water samples as small as 0.1 μL .

Analytical Method

High temperature conversion in continuous flow mode became available with the introduction of the glassy carbon reactors. The conversion is based on the reaction:

$$\text{H}_2\text{O} + \text{C} \rightarrow \text{H}_2 + \text{CO}^1$$

Water samples are entrained in a He carrier gas stream and passed through a glassy carbon reactor at 1400 °C. The tube-in-tube principle (glassy carbon tube in a ceramic tube, Figure 1) prevents oxidation of the glassy carbon by keeping the tube in an oxygen-free atmosphere, constantly flushed by a flow of He. Contamination or exchange of the O in the sample CO with O in the surrounding ceramic (Al_2O_3) tube is prevented because the reaction takes place inside the glassy carbon tube.

Isotope Ratio Analysis of Aqueous Solutions, e.g.

- water
- waste water
- beer
- wine
- fruit juice
- urine
- blood plasma
- ice
- atmospheric water


Figure 1: Schematic of the Thermo Scientific high temperature conversion elemental analyzer TC/EA.

The products of the reduction reaction, H₂ and CO, are separated on a gas chromatographic column and analyzed directly in a Thermo Scientific DELTA V Isotope Ratio Mass Spectrometer. After collecting the ion beams of *m/z* 2 and *m/z* 3 in the H₂ peak a fast switch of the magnetic field allows the subsequent collection of *m/z* 28 and *m/z* 30 in the CO peak (dual measurement mode). A narrow stainless steel insert is used at the top of the reactor to ensure both rapid and low memory transfer of the water into the hot zone. The system settings are given in Table 1 and 2.

It is essential that the injection of the water sample takes place immediately after the needle penetrates the septum to avoid isotopic fractionation resulting from evaporation in the needle tip. After the sample injection, the needle stays in the hot reactor for a certain period of time in order to remove any sample residue prior to the next injection.

Reactor temperature	1400 °C
GC temperature	100 °C
Carrier flow	100 ml/min
Autosampler type	Thermo Scientific AS 3000
Syringe size	0.5 µL

Table 1: High temperature conversion elemental analyzer (TC/EA) settings for H and O isotope ratio analysis.

Injection volume	0.1 µL
Plunger strokes	3
Rinses	0
Pre-injection dwell time	0 s
Post-injection dwell time	16 s
Solvent wash cycle	0

Table 2: AS 3000 Autosampler settings.

	$\delta^2\text{H}$ [‰] raw	Mean	S.D.	$\delta^2\text{H}$ [‰] after scaling	IAEA accepted value	$\delta^{18}\text{O}$ [‰] raw	Mean	S.D.	$\delta^{18}\text{O}$ [‰] after scaling	IAEA accepted value
GISP	-185.03					-24.23				
	-186.95					-24.67				
	-186.37					-24.75				
	-187.09					-24.75				
	-186.62	-186.76	0.32	-188.99	-189.5	-24.84	-24.75	0.07	-24.73	-24.5
SLAP	-420.53					-55.16				
	-422.13					-55.36				
	-422.35					-55.42				
	-424.45					-55.89				
	-422.67	-422.90	1.06	-428.00	-428.0	-55.52	-55.54	0.24	-55.50	-55.5
SMOW	-1.76					-0.32				
	-0.29					-0.04				
	0.31					-0.06				
	-0.02					-0.02				
	-0.11	-0.03	0.25	0.00	0.0	0.13	0.01	0.09	0.00	0.0
Sample A	-47.10					-7.44				
	-42.47					-6.99				
	-41.40					-6.82				
	-41.32					-6.66				
	-40.11	-41.32	0.96	-41.79	-	-6.76	-6.81	0.14	-6.81	-

Table 3: Mean and Standard Deviation (S.D.) of $\delta^2\text{H}$ and $\delta^{18}\text{O}$ values of five consecutive injections of international water standards GISP, SLAP, SMOW and a lab standard. The data is given in sequence of injection to demonstrate memory effects. The table shows raw data and data after scaling with SMOW and SLAP. The data was acquired during field installation. Strike-through data is not considered.

Results

Accuracy

Table 3 shows data from an installation of a TC/EA at a customer site. International water standards (IAEA, Vienna) plus a lab standard (Sample A) were injected in the given sequence. The water standards cover the full range of $\delta^2\text{H}$ and $\delta^{18}\text{O}$ values of natural water. The specifications of the TC/EA are < 2‰ for $\delta^2\text{H}$ and < 0.2‰ for $\delta^{18}\text{O}$. The typical precision achieved during routine analysis can be < 1‰ for $\delta^2\text{H}$ and < 0.1‰ for $\delta^{18}\text{O}$. All samples were analyzed in dual measurement mode.

The uncorrected raw data in Table 3 show that there is no need for any post-acquisition evaluation or correction. The direct conversion of the sample in the reactor, together with a highly precise, small and stable H₃⁺-factor of 6.2 ppm/nA provides very accurate delta values. This is a clear advantage over equilibration techniques. Common approaches for post-evaluation of raw delta values like the SMOW/SLAP scaling and drift corrections have been discussed by several authors.^{2, 3, 4}

Memory

In Table 3 only the first data point of each set of five injections shows the influence of preceding samples. Large differences between delta values of adjacent samples can corrupt the first analysis due to memory effects mainly caused by the syringe. The implementation of post-injection dwell time of the syringe in the reactor removes most of this memory. Washing cycles for the syringe do not help here because any solvent would contaminate the syringe with H- and O-isotopes.

Analysis Time

Table 4 shows a sequence of injections with samples of Antarctic precipitation with similar isotope ratio values. Mean and Standard Deviation (S.D.) were calculated from five consecutive injections per sample. The mean precision for $\delta^2\text{H}$ is $\ll 0.6\text{‰}$ which is a requirement for ice core analysis by paleoclimatologists. High resolution data of long ice cores demand high sample throughput.

Figure 2 and Figure 3 compare the mean value with the first injection of each set to see if the number of repetitions can be reduced without affecting the accuracy. All $\delta^{18}\text{O}$ isotope values of the first injection fall within the S.D. except for Sample H showing that already the first injection gives correct results. About half of the $\delta^2\text{H}$ values of the first injection fall outside the S.D. range requiring repetition.

Sample	Mean $\delta^2\text{H}$	S.D.	Mean $\delta^{18}\text{O}$	S.D.
Standard	-413.82	0.40	-39.48	0.07
Sample A	-429.78	0.43	-41.46	0.03
Sample B	-423.93	0.79	-40.68	0.09
Sample C	-422.30	0.55	-40.52	0.08
Sample D	-418.52	0.55	-40.02	0.06
Standard	-414.56	0.28	-39.58	0.11
Sample E	-402.63	0.48	-37.94	0.04
Sample F	-398.87	0.25	-37.54	0.08
Sample G	-392.22	0.88	-36.56	0.25
Sample H	-394.69	0.48	-36.83	0.13
Standard	-412.87	0.76	-39.15	0.13

Table 4: Mean and Standard Deviation (S.D.) of $\delta^2\text{H}$ and $\delta^{18}\text{O}$ values of five consecutive sample injections ($N = 5$) with similar isotope ratio values. No data was excluded as memory effects are negligible. Standard: $N = 3$.


Figure 2: $\delta^2\text{H}$ mean (sample $N = 5$, Standard $N = 3$) with Standard Deviation. Red circle data points represent the first injection of each set, all first injections fall within the S.D. range except for Sample H.


Figure 3: $\delta^{18}\text{O}$ mean (sample $N = 5$, standard $N = 3$) with Standard Deviation. Red circle data points represent the first injection of each set. Most first injections fall within the S.D. range.

In Short

Single Measurement $\delta^{18}\text{O}$

6 minutes analysis time,
no repetitions for screening required ($N = 1$),
10 injection per hour,
240 samples per day

Single Measurement $\delta^2\text{H}$

6-7 minutes analysis time,
2 repetitions (multiple injections, $N = 3$),
9 samples per hour,
216 samples per day
See Figure 4

Dual Measurement $\delta^2\text{H}$ and $\delta^{18}\text{O}$

5-6 minutes analysis time,
2 repetitions ($N = 3$),
9 injections per hour,
72 samples per day
(9 injections \times 24 h/3(triplicate) = 72)
See Figure 5


Figure 4: Single measurement $\delta^2\text{H}$ with multiple injections of the same sample to account for two repetitions. The analysis time is six minutes.


Figure 5: Dual measurement of $\delta^2\text{H}$ and $\delta^{18}\text{O}$ from a single injection with a fast magnet switch at 120 s. The analysis time can be reduced to less than five minutes.

Sample	Mean $\delta^2\text{H}$	S.D.	Mean $\delta^{18}\text{O}$	S.D.
Lab Standard	-53.69	0.20	-7.84	0.04
V-SMOW	0.00	0.34	0.00	0.08
GISP	-189.50	0.40	-24.80	0.05
Ethanol	-234.80	0.12	-24.18	0.08
Lab Standard	-53.76	0.39	-7.85	0.03

Table 5: Mean and Standard Deviation (S.D.) of $\delta^2\text{H}$ and $\delta^{18}\text{O}$ values of two standards and ethanol demonstrating the possibility to run organic liquids. N = 10 for standards, N = 5 for ethanol.

Other Liquid Samples

While isotope ratio analysis of water is by far the predominant application, other substrates have been tested, including purified blood plasma, urine, and ethanol (see Table 5), see also Calderone et al. for ethanol applications.⁵ The application lab of Thermo Fisher Scientific, Bremen, Germany, took part in an interlaboratory study initiated by BEVABS (Ispra), IRMM Geel, and the French bioanalytical company Eurofins. The object of the interlaboratory study was to evaluate the use of tetramethylurea (TMU) as a reference compound for TC/EA-IRMS analysis of wine and alcohols as defined in European Commission Regulation (EEC) N° 2676/90.

The final report of the TC/EA-IRMS interlaboratory study states the applicability and emerging potential of EA-IRMS for food authenticity testing.⁶ Even the analysis of organic substances like dodecane, methyl dodecanoate and methyl N-methylantranilate produced good results.

Conclusion

Since its introduction in 1997, the Thermo Scientific high temperature conversion elemental analyzer TC/EA has created a new field of application for the O- and H-isotope analysis world-wide. TC/EA-IRMS is since then used for the analysis of liquids.

The advantage of fast sub- μL analysis of liquid samples (only five minutes for two isotope ratios) by far outweighs the eventually lower precision.

Low investment with high sample throughput leads to low cost per analysis. In addition, an autosampler for solids with minor modifications on the reactor packing makes the TC/EA also a valuable system for the isotope ratio analysis of solid samples.

Since 2003, the oven design and reactor technology of the TC/EA has been integrated in the Thermo Scientific Flash 2000 HT elemental analyzer for isotope ratio MS, so that Dumas combustion for the analysis of N and C isotopes is combined with high temperature conversion for analysis of O and H isotopes in a single system.

The Flash 2000 HT elemental analyzer for isotope ratio MS enables the user to switch from N-C isotope ratio analysis to H-O isotope ratio analysis within just a few minutes. The performance of water analyses in both elemental analyzers the TC/EA and the Flash 2000 HT is identical.

References and Related Literature

- 1 Koziat J.J., Isotope Ratio Mass Spectrometric Method for the Online Determination of Oxygen-18 in Organic Matter, *J. Mass Spectrom* **1997**, *32*, 103-108.
- 2 Coplen T., Normalization of oxygen and hydrogen isotope data, *Chem. Geol.* **1988**, *72*, 293-297.
- 3 Nelson S.T., A simple, practical methodology for routine VSMOW/SLAP normalization of water samples analyzed by continuous flow methods, *RCM* **14**, **2000**, 1044-1046.
- 4 Werner R. & Brand W.A., Referencing strategies and techniques in stable isotope ratio analysis, *RCM* **15**, **2001**, 501-519.
- 5 Calderone G., Reniero F. & Guillou C., $^{18}\text{O}/^{16}\text{O}$ measurements on ethanol, *RCM* **20**, **2006**, 937-940.
- 6 Bréas O., Thomas F., Zeleny R., Calderone G., Jamin E. & Guillou C., Performance evaluation of elemental analysis/isotope ratio mass spectrometry methods for the determination of the D/H ratio in tetramethylurea and other compounds – results of a laboratory inter-comparison, *RCM* **21**, **2007**, 1555-1560.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa-Other
+27 11 570 1840

Australia
+61 3 9757 4300

Austria
+43 1 333 50 34 0

Belgium
+32 53 73 42 41

Canada
+1 800 530 8447

China
+86 10 8419 3588

Denmark
+45 70 23 62 60

Europe-Other
+43 1 333 50 34 0

Finland/Norway/Sweden
+46 8 556 468 00

France
+33 1 60 92 48 00

Germany
+49 6103 408 1014

India
+91 22 6742 9434

Italy
+39 02 950 591

Japan
+81 45 453 9100

Latin America
+1 561 688 8700

Middle East
+43 1 333 50 34 0

Netherlands
+31 76 579 55 55

New Zealand
+64 9 980 6700

Russia/CIS
+43 1 333 50 34 0

South Africa
+27 11 570 1840

Spain
+34 914 845 965

Switzerland
+41 61 716 77 00

UK
+44 1442 233555

USA
+1 800 532 4752

Thermo Fisher Scientific (Bremen)
GmbH Management System
Registered to ISO 9001:2008

AN30180_E 11/10C

www.thermoscientific.com

©2010 Thermo Fisher Scientific Inc. All rights reserved. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.