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Operating Conditions	Water phase	Organic Phase
	HNO <sub>3</sub> /HCl	DMSO
<b>Sample Introduction System</b>		
Nebulizer	Quartz	PFA
Injector	Quartz 2.5mm ID	Quartz 1.0mm ID
Spray Chamber	Quartz, cyclonic	Quartz, cyclonic
Interface	Pt cone	Pt cone
<b>Parameter</b>		
Plasma mode	KED	KED
RF forward power (w)	1550	1350
Sampling depth (mm)	5	7.5
Nebulizer gas Flow (L/min)	1.03	0.688
Spray Chamber Temperature(°C)	2.73	20
He cell gas flow (mL/min)	4.2	4
Addition gas flow (mL/min)	none	O <sub>2</sub> , 30

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# Online Multi-elemental Monitoring of Environmental Atmospheric Gases with a Gas Exchange Device Coupled to the High Sensitivity Thermo Scientific iCAP Qs ICP-MS

Tomoko Vincent<sup>1</sup>, Kohei Nishiguchi<sup>2</sup>, Keisuke Utani<sup>2</sup>, Shona McSheehy Ducos<sup>1</sup>, Lothar Rottmann<sup>1</sup>

<sup>1</sup>Thermo Fisher Scientific, Hanna-Kunath Strasse 11, 28199 Bremen, Germany

<sup>2</sup>J-SCIENCE LAB, 3-1 Hiuchigata, Kamitoba, Minami-ku, 601-8144 Kyoto, Japan

## Key Words

Environmental gas monitoring, KED, Nuclear, Air, Pollution

## Goal

To demonstrate real time multi-elemental analysis of environmental contaminants in air in using a gas exchange device coupled to the Thermo Scientific™ iCAP™ Qs ICP-MS.

## Introduction

Environmental pollution is an immediate and growing concern. Environmental monitoring has become a key mechanism in determining how industrial activities and accidental contamination impact our water supplies, ecosystems and the air that we breathe. Strategic policies such as Clean Air for Europe (CAFE) and the United States Clean Air Act (CAA) for example, provide a framework of standards and objectives for the control of persistent pollutants that can damage health and the environment. Daughter directives within the CAFE program call for monitoring of nickel, arsenic, cadmium, mercury and lead. Of particular concern are fine particles with a diameter of less than 2.5  $\mu\text{m}$  (PM<sub>2.5</sub>), which can penetrate deep into human bronchial tubes causing asthma and bronchitis. In China, where coal burning is still the biggest energy source and private car use is rapidly increasing, an estimated 70% of cities have PM<sub>2.5</sub> levels above<sup>1</sup> the daily 75  $\mu\text{g}$  PM<sub>2.5</sub>/m<sup>3</sup> Chinese National Air Quality Standard limit<sup>2</sup>.

The Fukushima Daiichi nuclear power plant incident in 2011 demonstrated an immediate need for the monitoring of specific analytes in ambient air<sup>3</sup> in order to determine the distribution of radioactive materials in the environment, support remediation strategies and assess any possible subsequent threat to human health.

Current approaches for measuring radioactive contamination in the environment require the collection and preparation of soil and water samples and impingers or filtering for air sampling. These methods, however, cannot provide the real time results necessary for a meaningful assessment of current environmental conditions.

This application note evaluates the use of a gas exchange device (GED) coupled to the iCAP Qs ICP-MS for the direct analysis of atmospheric air (Figure 1). The GED overcomes problems in the direct analysis of air by ICP based techniques by exchanging atmospheric gases with argon which is compatible with the ICP ion source. Direct air sampling using the GED and on-line elemental analysis by high sensitivity ICP-MS therefore provides immediate information on airborne contamination.

The fast elemental scan speeds afforded by ICP-MS allow for the measurement of single particle events (SPE) that are seen as pulses in the signal intensity as they are processed by the ICP ion source. As air entrained elements are in particulate form, the ICP-MS based analysis of SPE can determine the number and even size of particles, therefore providing valuable information on elemental transport in the environment.

## Gas Exchange Device

Figure 1. Gas Exchange Device coupled to the Thermo Scientific iCAP Qs

The GED (J-SCIENCE LAB Co. Ltd., Kyoto Japan) consists of 2 inner tubes of 0.07  $\mu\text{m}$  porous glass that act as a membrane and an outer body made of (PYREX™) glass (Figure 2). The gas sample is introduced into the central inner tube of the GED device and argon is introduced into the outer tube at 2 L/min using the GED internal mass flow controller.

As the particulate supporting gas travels along the inner tube, atmospheric gases diffuse out of the inner tube (indicated by the blue arrows in Figure 2) and argon diffuses from the outer tube into the inner tube (indicated by the yellow arrows). Ultimately, before the sample enters the plasma, the atmospheric gases have been replaced by argon while particulate matter from the sample remains.


Figure 2. Principle of Gas Exchange Device

## Instrument Configuration

A Thermo Scientific iCAP Qs ICP-MS was used for all measurements. The sample introduction kit consisted of a demountable quartz torch with a 2.0 mm ID sapphire injector. The spray chamber was removed and the exit line from the GED was attached directly to the injector using a PFA adapter. The instrument was operated in a single QCell mode, using kinetic energy discrimination (KED) with pure He as the collision gas. The use of pure He KED eliminates interferences, ensuring low backgrounds and maintains high sensitivity for the analytes measured.

Table 1. Instrument operating parameters

Parameter	Value
Sample gas	0.68 L/min
Air sampling rate	0.235 L/min
Auxiliary gas	0.8 L/min
Cool gas	14 L/min
Forward power	1550 W
Collision cell gas	He at 4.5 mL/min
KED barrier	3 V

## Tuning, Sampling and Data Acquisition

The combined GED-ICP-MS system was optimized using a gas mixture of volatile  $^{52}\text{Cr}$ ,  $^{95}\text{Mo}$ , and  $^{184}\text{W}$  species generated by an Element Standards Gas Generator (ESGG) (J-SCIENCE LAB). The ESGG produces a continuous stream of elemental vapor in an argon gas flow that is introduced into the GED.

Ambient air was sampled via a TYGON® R3603 tube (i.d. 3.2mm, o.d. 6.4mm), approximately 15 m long. One end of the tube was connected directly to the GED; the other end was fed through an outside window, 10 m above ground. The argon sample gas flow of 0.68 L/min entrained the outside air through the tube and into the GED at a flow rate of 0.235 L/min. Sampling conditions are shown in Table 2.

Table 2. Sampling Conditions

Parameter	Value
Location	Bremen, Germany
Date, Time	18 Dec 2012, 9:30 am to 13:30 pm
Weather Conditions	Overcast, 6 °C, Wind: ENE 14 km/h

To determine optimal dwell times for the analysis of atmospheric particles, time resolved analyses of the target isotopes ( $^{66}\text{Zn}$ ,  $^{118}\text{Sn}$ ,  $^{121}\text{Sb}$ ,  $^{141}\text{Pr}$ ,  $^{205}\text{Tl}$ ,  $^{206, 207, 208}\text{Pb}$ ,  $^{209}\text{Bi}$  and  $^{238}\text{U}$ ) were made using dwell times of 0.5, 1, 10 and 100 ms over a period of 300 s. In a subsequent 4 hour measurement of the outside air, a 10 ms dwell time was used for all isotopes.

## Result and discussion

Ambient air was self-aspirated into the GED-ICP-MS and the target isotopes were analyzed at different dwell times as described. For analyses that intend to evaluate particle information, dwell times should be carefully optimized to favor the analysis of full SPEs. Dwell times that are too short, will measure sections of the SPE and dwell times that are too long will tend to measure multiple particle events in one scan.

Figure 3 (100 ms dwell) and Figure 4 (10 ms dwell) show time resolved scans for both an N<sub>2</sub> gas blank and an air sample. When comparing Figures 3 and 4, it is clear that

10 ms dwell time (higher scan frequency) provides improved discrimination between particle events and that individual events have a higher signal to noise ratio. At 100 ms dwell time, multiple events are captured together and averaged so that particle event information such as the number and intensity of events is lost.

Even shorter dwell times should further improve time resolution and dwell times of 0.5 and 1 ms were also evaluated. It was observed however, that shorter dwell times captured partial particle events and did not provide additional information compared to the 10 ms dwell time scans.


Figure 3. Time resolved scans for ambient air and an N<sub>2</sub> gas blank using a 100 ms dwell time. Please note the different (vertical) intensity scales on the two scans


Figure 4. Time resolved scans for outside air and an N<sub>2</sub> gas blank using a 10 ms dwell time. Please note the different (vertical) intensity scales on the two scans

A 10 ms dwell time was chosen for a longer term measurement of over 4 hours. Figure 5 shows a 30 minute section of the 4 hour monitoring period. Particle events for Zn, Pb and Bi are clearly shown.


Figure 5. 30 minute section of the 4 hour air monitoring period

The environmental source of the particles observed in these analyses is practically impossible to determine as the particle trajectory travel distance is unknown. In this application, lead particles can be seen at regular intervals throughout the 4 hour period. The cycling and transport of lead in particulate matter is well documented, therefore lead particles are likely to occur in industrial and urban areas. For a more accurate evaluation of particle sources in air, further sampling points and baseline data regarding the geochemical cycling, potential sources and detailed environmental conditions at the time of sampling are required.

## Conclusion

The Thermo Scientific iCAP Qs ICP-MS has demonstrated the high sensitivity and freedom from background interferences required for the measurement of particles in ambient air. The iCAP Qs ICP-MS has flexible scan times with a data buffer system capable of handling the large amounts of data acquired during long term *in-situ* analyses. QCell technology eliminates interferences whilst ensuring high sensitivity even for the lower mass isotopes thanks to its high transmission efficiency. A single mode KED approach could thus be employed for a number of isotopes across the mass range.

The use of a Gas Exchange Device coupled to the iCAP Qs ICP-MS is not limited to monitoring environmental air. Additional applications include monitoring a large range of gas samples such as specialty gases, semiconductor gases in the workplace, tobacco smoke and the presence of radionuclides.

## Acknowledgement:

The gas exchange device and gas reaction system were provided by J-SCIENCE Lab., Japan

1. Xinhuanet News [http://news.xinhuanet.com/english/china/2012-12/05/c\\_132021534.htm](http://news.xinhuanet.com/english/china/2012-12/05/c_132021534.htm)
2. Beijing Environmental Protection Bureau (BJEPB) (2012) <http://www.bjepb.gov.cn>
3. Carlsbad Environmental Monitoring and Research Center 2011 Annual Report [www.cemrc.org/wp-content/uploads/2011/09/2011-Annual-Report.pdf](http://www.cemrc.org/wp-content/uploads/2011/09/2011-Annual-Report.pdf)

[www.thermoscientific.com](http://www.thermoscientific.com)

©2013 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Thermo Fisher Scientific, (Bremen) GmbH  
Management System Registered is  
ISO 9001:2008 Certified.

<b>Africa-Other</b> +27 11 570 1840	<b>Europe-Other</b> +43 1 333 50 34 0	<b>Japan</b> +81 45 453 9100	<b>Spain</b> +34 914 845 965
<b>Australia</b> +61 3 9757 4300	<b>Finland/Norway/Sweden</b> +46 8 556 468 00	<b>Latin America</b> +1 561 688 8700	<b>Switzerland</b> +41 61 716 77 00
<b>Austria</b> +43 1 333 50 34 0	<b>France</b> +33 1 60 92 48 00	<b>Middle East</b> +43 1 333 50 34 0	<b>UK</b> +44 1442 233555
<b>Belgium</b> +32 53 73 42 41	<b>Germany</b> +49 6103 408 1014	<b>Netherlands</b> +31 76 579 55 55	<b>USA</b> +1 800 532 4752
<b>Canada</b> +1 800 530 8447	<b>India</b> +91 22 6742 9434	<b>New Zealand</b> +64 9 980 6700	
<b>China</b> +86 10 8419 3588	<b>Italy</b> +39 02 950 591	<b>Russia/CIS</b> +43 1 333 50 34 0	
<b>Denmark</b> +45 70 23 62 60		<b>South Africa</b> +27 11 570 1840	

# Analysis of plating baths using the Thermo Scientific iCAP 7400 ICP-OES

Patrícia Coelho, Applications Chemist, Thermo Fisher Scientific, Cambridge, UK  
Grégory Lecornet, Applications Chemist, Thermo Fisher Scientific, Courtaboeuf, France

## Key Words

Method of standard addition, plating baths, sulphuric acid levels

## Goal

This application note describes the performance of the Thermo Scientific™ iCAP™ 7400 ICP-OES Duo analysing different elements in different types of plating baths. The duo viewing offers optimal method conditions using axial view for traces and radial view for major elements, and the pre-loaded template allows for rapid and simple method development.

## Introduction

The plating of metal is an ancient technique, used for hundreds of years. It can be defined as the act of covering the surface of objects by depositing a metal on a conductive surface. This is typically carried out by immersing the object in a solution in which the metal ions are moved by an electric field. The applications are rather vast, ranging from the merely decorative, to the enhancement of the physical properties of the material being covered (i.e. to prevent corrosion, reduce friction, alter conductivity, or to improve characteristics such as hardness or durability). It is widely used in the production of jewellery pieces to achieve a silver or gold finish. Due to the ability to cover objects as small as an atom this technique has a high potential to be used in the nanotechnology area.

Manufacturers frequently establish optimum specifications in order to guarantee maximum solution efficiency and uniformity of the plating solution. Plating solutions should be analysed regularly in order to maintain the recommended formulations and to prevent problems related to improper levels of bath constituents and contaminants. A current problem faced by the plating industry are gradual and continuous changes in pH and metal or cyanide content, leading to a significant decrease in efficiency. Some of the most common plating solutions and their main applications are presented in table 1.

Table 1. Plating solutions commonly used and their main applications

Plating Baths	Main applications
Alloy	Hardness improvement
Cadmium	Corrosion Resistance
Chrome	Decorative and industrial
Composite	Physical properties improvement
Golden	Jewellery manufacturing and electronics
Nickel	Decorative and corrosion resistance
Rhodium	Jewellery manufacturing
Silver	Jewellery manufacturing and electronics
Tin	Food processing and electronics
Zinc	Corrosion resistance

There are several different methods that can be employed for the quantitative analysis of plating solutions; these are classified as volumetric, gravimetric or instrumental. Both volumetric and gravimetric methods are simple, accurate, and rapid and can be performed with common laboratory equipment. Nevertheless they only rely on chemical reactions instead of measuring the physical properties associated with the composition of the substance. Additionally, instrumental methods are far quicker and allow for the automatization of the analysis, leading to less mathematical errors and higher reproducibility.

A common instrumental technique is spectroscopic analysis, in particular inductively coupled plasma-optical emission spectroscopy (ICP-OES), which is used in the analysis of major components and trace contaminants in plating solutions. When using this technique there are some aspects that need to be taken into consideration, such as physical interferences (such as viscosity or surface tension) and chemical interferences. They can be overcome easily by sample dilution and accurate matrix matching.

The analysis of sulphuric acid ( $H_2SO_4$ ) levels in plating baths is commonly performed by electroanalytical methods, specifically potentiometry. This is a simple and relatively low cost method but the sensitivity is limited at low concentrations, not allowing for accurate measurements in highly dilute solutions and it can also encounter a number of interferences. Manufacturers that use plating in their processes, are increasingly using ICP-OES instrumentation to quantify  $H_2SO_4$  levels. This method allows for a rapid, sensitive and interference free measurement of sulphur, directly proportional to  $H_2SO_4$ , resulting in high accuracy even at low levels.

## Instrumentation

The Thermo Scientific iCAP 7400 ICP-OES Duo was used for the analysis of a range of plating baths. This is a compact duo view ICP-OES instrument based on the innovative technologies of the iCAP 7000 Series ICP-OES spectrometers. The instrument achieves powerful analyte detection and provides a highly cost effective solution for routine analysis of liquids in laboratories with standard sample throughput requirements. The Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ (ISDS) incorporates several pre-loaded templates (see Figure 1) for common methods, simplifying normal method development and providing an option of immediate analysis.


Figure 1. Plating Baths template selection

## Sample and standard preparation

Elements analysed in each plating bath are presented in table 2. For the Chromium and Chromium & Aluminium baths the quantification of sulphur was performed in order to quantify the amount of  $H_2SO_4$  in the samples.

Table 2. Elements analysed in each plating bath

Type of plating bath	Elements analysed
Chromium and ammonium fluoride (Cr & $NH_4F$ )	Al, Cr, Zr
Nickel (Ni)	B, Ni
Chromium (Cr)	Fe, S
Chromium (Cr)	S
Chromium and aluminium (Cr & Al)	S

Quantification of the samples was made by Method of Standard Addition (MSA) to avoid any possible matrix effect. Standard addition calibration was prepared by spiking samples with traceable 1000mg/L aqueous, single element standards, to the concentrations listed in table 3. All solutions were made to 50 mL with ultra pure de-ionized water.

Table 3. Standard addition concentrations (mg/L)

Type of Plating Bath	Element	Blank	Standard 1	Standard 2	Standard 3
Cr & NH <sub>4</sub> F	Al	0	1	2	-
	Cr	0	20	40	-
	Zr	0	50	100	-
Ni	B	0	5	10	-
	Ni	0	50	100	-
Cr	Fe	0	50	100	-
	S	0	50	100	-
Cr	S	0	50	100	150
Cr & Al	S	0	25	50	75

## Method Development

A LabBook was created in Qtegra ISDS which contained the method parameters and standard concentrations as listed in this note. A standard sample introduction kit was used for the analysis. The instrument was calibrated and the samples analyzed in a single run. The method parameters are shown below in table 4.

Table 4. Method parameters

Parameter	Setting
Pump tubing	Sample Tygon® 1.016 mm Drain Tygon® 1.524 mm
Pump speed	50 rpm
Nebulizer	Glass concentric
Nebulizer gas flow	0.5 L/min
Spray chamber	Glass cyclonic
Center tube	2 mm
RF Power	1150 W
Coolant gas flow	12 L/min
Auxiliary gas flow	0.5 L/min
Exposure times	Low 10 sec High 5 sec

## Results

The results obtained in the analysis of the different bath samples are shown in table 5. It is assumed that all of the sulphur present in the sample is in the form of H<sub>2</sub>SO<sub>4</sub>, therefore the results obtained for sulphur is multiplied by 3.06 [ $M(H_2SO_4)/M(S)=3.06$ ], in order to calculate the concentration of H<sub>2</sub>SO<sub>4</sub>. Results obtained for all the elements were within the expected range.

Table 5. Results of the analysis of the different plating baths. All concentrations are in mg/L.

Type of plating bath	Elements and wavelengths (nm)	Concentration found	Established concentration/range	Dilution factor*
Cr & NH <sub>4</sub> F	Al 167.079	2.2	<10	100
	Cr 205.560	456.0	500	
	Cr 206.550	454.9		
	Cr 267.716	446.7		
	Zr 274.256	1504	1700	
	Zr 327.305	1486		
	Zr 339.198	1567		
	Zr 343.823	1689		
Ni	B 208.959	6380	<7000	1000
	B 249.678	6270		
	B 249.773	6150		
	Ni 221.647	95470	80000	
	Ni 231.604	94770		
Cr	Fe 259.837	2892	<5000	100
	Fe 259.940	2907		
	Fe 371.994	2823		
	S 180.731 (H <sub>2</sub> SO <sub>4</sub> )	1021.6 (3126)	2500 – 3000 (H <sub>2</sub> SO <sub>4</sub> )	
	S 182.034 (H <sub>2</sub> SO <sub>4</sub> )	1035.0 (3167)		
	S 182.624 (H <sub>2</sub> SO <sub>4</sub> )	983.7 (3010)		
Cr	S 180.731 (H <sub>2</sub> SO <sub>4</sub> )	85.6 (2020.3)	2250 – 2500 (H <sub>2</sub> SO <sub>4</sub> )	10
	S 182.034 (H <sub>2</sub> SO <sub>4</sub> )	79.4 (2429.6)		
	S 182.624 (H <sub>2</sub> SO <sub>4</sub> )	74.2 (2271.7)		
Cr & Al	S 180.731 (H <sub>2</sub> SO <sub>4</sub> )	5.7 (35.1)	<150 (H <sub>2</sub> SO <sub>4</sub> )	2
	S 182.034 (H <sub>2</sub> SO <sub>4</sub> )	7.5 (45.9)		
	S 182.624 (H <sub>2</sub> SO <sub>4</sub> )	7.5 (45.9)		

\* All dilutions were prepared using ultra pure water.

Concentration ranges were established by the manufacturers taking into account the concentrations usually found in these types of plating baths to obtain maximum efficiency. This gives the plating operators the status of the bath simplifying the assessment of how the bath composition evolves overtime.

## Conclusion

The analysis of plating baths, specifically those that require the analysis of sulphuric acid, is rapid and highly sensitive when using the Thermo Scientific iCAP 7400 ICP-OES in conjunction with the Qtegra ISDS. This enables the accurate quantification of metals and sulphuric acid content, while the powerful, easy to use iCAP 7400 ICP-OES allows both experienced and inexperienced users alike to vastly reduce the method development time required for these sample types, resulting in cost effective analyses.

[www.thermoscientific.com](http://www.thermoscientific.com)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. Tygon is a registered trademark of Saint-Gobain. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Thermo Electron Manufacturing Ltd  
(Cambridge & Courtaboeuf)  
is ISO Certified.

**Africa-Other** +27 11 570 1840  
**Australia** +61 3 9757 4300  
**Austria** +43 1 333 50 34 0  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** +86 10 8419 3588  
**Denmark** +45 70 23 62 60

**Europe-Other** +43 1 333 50 34 0  
**Finland/Norway/Sweden**  
+46 8 556 468 00  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9434  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Russia/CIS** +43 1 333 50 34 0  
**South Africa** +27 11 570 1840

**Spain** +34 914 845 965  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752

AN43150\_E 01/13C

**Thermo**  
SCIENTIFIC  
Part of Thermo Fisher Scientific